


Betatech[®]
Hernia Repair Solutions





Betatech[®]

Hernia Repair Solutions

	VENTRAL	UMBILICAL	INGUINAL
Polymesh Polypropylene -Ultralight -Light -Standard -Blue	□	□	□
Polymesh Hernia System			△
Polymesh 3D			△
Polymesh Plug&Patch			△
Polymesh Composite	□	□	□
Polymesh Dual	○	○	
Polymesh Inova	○	○	

Polymesh

Hernia Repair Solutions

- Flat Mash
- △ Mesh Device
- Tissue Seperating Mesh



MULTI-USE

Polymesh

Polypropylene

Knitted Polypropylene
Monofilament Construction
Superior Biocompatibility and
Healing Response
Various weight and size range




Betatech[®]

Polymesh Polypropylene

Polypropylene mesh for the repair of abdominal wall fascial defects

	Size	Product Code				Units per Box
		Ultralight (28g/m ²)	Light (50g/m ²)	Standard (85g/m ²)	Blue (50g/m ²)	
	6x11cm	PMUL611	PM611	PMS611	PMB611	5 Per Box
	7,5x15cm	PMUL7515	PM7515	PMS7515	PMB7515	5 Per Box
	10x15cm	PMUL1015	PM1015	PMS1015	PMB1015	5 Per Box
	15x15cm	PMUL1515	PM1515	PMS1515	PMB1515	5 Per Box
	30x30cm	PMUL3030	PM3030	PMS3030	PMB3030	5 Per Box
	4,5x10cm	PMUL4510	PM4510	PMS4510	PMB4510	5 Per Box
	5x11cm	PMUL511	PM511	PMS511	PMB511	5 Per Box
	5x12cm	PMUL512	PM512	PMS512	PMB512	5 Per Box
	6x13cm	PMUL613	PM613	PMS613	PMB613	5 Per Box
	4,5x10cm	PMULE4510	PME4510	PMSE4510	PMBE4510	5 Per Box
	5x11cm	PMULE511	PME511	PMSE511	PMBE511	5 Per Box
	5x12cm	PMULE512	PME512	PMSE512	PMBE512	5 Per Box
	6x13cm	PMULE613	PME613	PMSE613	PMBE613	5 Per Box
	Diameter 5cm	PMULO5	PMO5	PMSO5	PMBO5	5 Per Box
	Diameter 7cm	PMULO7	PMO7	PMSO7	PMBO7	5 Per Box
	Diameter 5cm	PMULE5	PME5	PMSE5	PMBE5	5 Per Box
	Diameter 7cm	PMULE7	PME7	PMSE7	PMBE7	5 Per Box
	Diameter 5cm	PMULOE5	PMOE5	PMSOE5	PMBOE5	5 Per Box
	Diameter 7cm	PMULOE7	PMOE7	PMSOE7	PMBOE7	5 Per Box

Hernia Repair Solutions





INGUINAL

Polymesh

Plug&Patch

Mesh plug designed for use in a tension-free preperitoneal open inguinal hernia repair technique.







Betatech[®]

Polymesh Plug&Patch

Polymesh Plug&Patch is designed with pleated edges that conform readily to defects. The variety of sizes and models accommodates all groin hernias.



	Product Code	Size / Model	Units per Box
	PMP2	Plug H:2cm	1 Per Box
	PMP3	Plug H:3cm	1 Per Box
	PMP4	Plug H:4cm	1 Per Box
	PMP5	Plug H:5cm	1 Per Box
	<hr/>		
	PLP01	Plug&Patch H:3,8cm, 6x12cm	1 Per Box
	PLP02	Plug&Patch H:4,2cm, 6x12cm	1 Per Box
<hr/>			
	PLPC01	Plug&Patch (With Part.Absorbable Patch)	1 Per Box
	PLPC02	Plug&Patch (With Part.Absorbable Patch)	1 Per Box





INGUINAL

Polymesh

Hernia System

Pre-shaped, three dimensional mesh constructed of an onlay patch connected by a cylinder to a circular underlay patch.



Betatech[®]


Polymesh Hernia System

Polymesh Hernia System has a bilayer design proven to combine anterior and posterior repair

Product Code	Size / Model	Units per Box
PPHSM	Diameter 8cm	1 Per Box
PPHSL	Diameter 10cm	1 Per Box

Hernia Repair Solutions





INGUINAL

Polymesh 3D


Anatomically
Shaped Polypropylene Mesh
Constructed of knitted
polypropylene monofilaments.






Betatech[®]

Polymesh 3D

Eliminates need for fixation, which saves time and cost.
Reduced pain for the patients who received Polymesh 3D compared with received flat mesh.



Product Code	Size / Model	Units per Box
PM3DM-L	3D Medium Left	1 Per Box
PM3DL-L	3D Large Left	1 Per Box
PM3DXL-L	3D X-Large Left	1 Per Box
PM3DM-R	3D Medium Right	1 Per Box
PM3DL-R	3D Large Right	1 Per Box
PM3DXL-R	3D X-Large Right	1 Per Box
PML3DM-L	3D Light Medium Left	1 Per Box
PML3DL-L	3D Light Large Left	1 Per Box
PML3DXL-L	3D Light X-Large Left	1 Per Box
PML3DM-R	3D Light Medium Right	1 Per Box
PML3DL-R	3D Light Large Right	1 Per Box
PML3DXL-R	3D Light X-Large Right	1 Per Box





INGUINAL

Polymesh

Composite


Partially Absorbable
Composite Mesh
manufactured from
absorbable PGA-PCL
and non-absorbable
polypropylene monofilament.



Betatech[®]

Polymesh Composite

Designed for patient comfort, comfortable healing. Absorbable in 90-120 days. Low-weight polypropylene PGA-PCL mesh with pore size of 2 to 4 mm provides less foreign body reaction and low postoperative adhesion.



Product Code	Size / Model	Units per Box
SCM0611	6x11cm	1 Per Box
SCM0614	6x14cm	1 Per Box
SCM7515	7,5x15cm	1 Per Box
SCM1012	10x12cm	1 Per Box
SCM1015	10x15cm	1 Per Box
SCM1515	15x15cm	1 Per Box
SCM1520	15x20cm	1 Per Box
SCM2020	20x20cm	1 Per Box
SCM2030	20x30cm	1 Per Box
SCM3030	30x30cm	1 Per Box



Before Partial Absorption



During Partial Absorption

Hernia Repair Solutions



VENTRAL

Polymesh

Dual

Dual Side Surgical Mesh
Silicone Coated Nonwoven
Polypropylene



Betatech[®]

Polymesh Dual

Constructed of double sides with different properties for implantation in intraperitoneal site. Visceral side is made of non-adherent silicone which prevents adhesion formation and creates a pseudo-peritoneum. Parietal side is made of nonwoven polypropylene for an easier tissue fixation and colonization.



Product Code	Size / Model	Units per Box
SDM1022	10x15/22cm (barrel shape)	1 Per Box
SDM1015	10x15cm	1 Per Box
SDMR12	Diameter 12cm	1 Per Box
SDMR15	Diameter 15cm	1 Per Box
SDMR18	Diameter 18cm	1 Per Box
SDM1515	15x15cm	1 Per Box
SDM1520	15x20cm	1 Per Box
SDM1530	15x30cm	1 Per Box
SDM1525	15x25cm	1 Per Box
SDM2020	20x20cm	1 Per Box
SDM2030	20x30cm	1 Per Box
SDM2535	25x35cm	1 Per Box
SDM2634	26x34cm	1 Per Box
SDM3030	30x30cm	1 Per Box



VENTRAL

Polymesh

Inova


Partially Absorbable Composite Mesh manufactured from absorbable PGA-PCL film coated non-absorbable polypropylene monofilament.




Betatech[®]

Polymesh Inova

Designed for open and laparoscopic hernia repair when a tissue separating layer is needed. POLYMESH INOVA combines Betatech's Polymesh Polypropylene with a tissue separating film layer of PGA-PCL absorbed within 90-120 days



Product Code	Size / Model	Units per Box
ISCM1022	10x15/22cm (barrel shape)	1 Per Box
ISCM1015	10x15cm	1 Per Box
ISCMR12	Diameter 12cm	1 Per Box
ISCMR15	Diameter 15cm	1 Per Box
ISCMR18	Diameter 18cm	1 Per Box
ISCM1515	15x15cm	1 Per Box
ISCM1520	15x20cm	1 Per Box
ISCM2030	20x30cm	1 Per Box
ISCM2634	26x34cm	1 Per Box
ISCM3030	30x30cm	1 Per Box





Hernia Repair Solutions

For additional ordering and size information, please visit www.betatechmedical.com,
contact your Medilution Healthcare representative or call +92 21 3483 2395



Ikitelli Org. San. Blg. Ataturk Oto San. Sit. 22. Sk.
Unal Is Merkezi No: 9, Basaksehir / Istanbul, TURKEY
P:+90 212 6710683 F: +90 212 6710684
info@betatechmedical.com
betatechmedical.com

MEDILUTION
HEALTH CARE



A-322 First Floor, Gulshan-e-Iqbal
Block 2, KDA Scheme No. 24,
Karachi.
P:+92 21 3483 2395
info@medilutionhc.com
www.medilutionhc.com